

Advanced technologies from



As Featured by Gilberto Debelian, DMD, PhD
in his Master Clinical Workshops and Lectures


ENDODONTIC FILE FORCEPS

- Reduces incidence of joint fatigue when using K-files
- Quick connect/disconnect
- Greater ergonomics
- Two distinct models to choose from:

PCF-N-90AHF/L




PCF-N-75CHF/L



The sharp prongs may impale the parachute holes facilitating 270°s of positioning in the X, Y and Z axis for easy approach to all canals.

As file forceps they are also ideal for use with mirrors or microscopes. For obturation they are ideal for use with Gutta-Percha and resin based carriers.



STEIGLITZ FORCEPS



North/South angulations address the bucco-lingual chamber orientation of anteriors and bicuspids.

Available in 45°, 75° and 90°

- With carbide inserts (extra grip)
- With diamond dust (extra thin)



East/West angulations address the mesio-distal chamber orientation of molars.

Available in 45°, 75° and 90°

- With carbide inserts (extra grip)
- With diamond dust (extra thin)