

THE
ULTIMATE
SHINE
IN LESS TIME



WHITE
JAZZ[™]
POLISHERS

SSWHITE® JAZZ™ POLISHERS



PROBLEM: What polishing system do I use on my many different composite materials that will achieve the highest shine possible without requiring multiple steps?

SOLUTION: jazz™ offers a polishing system that works with all composite materials and achieves the ultimate shine in a single step.

When only the highest shine will do:

JAZZ
SUPREME | COMPOSITE
1-STEP
UNIVERSAL
POLISHING SYSTEM

Breakthrough technology produces the highest shine possible.

Single-step system saves time.

Compatible with all direct aesthetic materials which eliminates the need for multiple polishing systems.

Re-usable for multi-patient-use.



Disc
Item #JAZZ



Cup
Item #JAZZ



Rock Edge
Item #JAZZ



Small Flange
Item #JAZZ

FOR COMPOSITES

PROBLEM: Choosing a composite polishing system is confusing since there seems to be a separate system for each material.

SOLUTION: Jazz™ offers a polishing system that is compatible with all composite materials.

When you need 1 polishing system for all composite materials:

JAZZ
C2S | COMPOSITE
2-STEP
POLISHING
SYSTEM

Diamond impregnated to achieve a high luster on all direct esthetic materials.

Compatible with all direct esthetic materials which eliminates the need for multiple polishing systems.

Re-usable for multi-patient-use.

1. Use the medium grit polisher for reducing and smoothing.



Medium Flute
Item #1002



Medium Cap
Item #1003



Medium Knit Edge
Item #1007



Medium Flat Flare
Item #1001

2. Finish with the fine grit polisher to achieve a high luster.



Fine Flute
Item #1005



Fine Cap
Item #1006



Fine Knit Edge
Item #1017



Fine Flat Flare
Item #1004

PROBLEM: How do I avoid the time and expense of cleaning and autoclaving composite polishers? Is there a system available that achieves a natural shine with the convenience of disposable instruments?

SOLUTION: Jazz™ offers a single-patient-use composite polishing system that creates a natural shine.

When you need to save time without compromising quality:

JAZZ
C1S | COMPOSITE
1-STEP
SINGLE-PATIENT
POLISHING SYSTEM

Diamond impregnated to achieve a high gloss on all direct esthetic materials.

Compatible with all direct esthetic materials which eliminates the need for multiple polishing systems.

Single-patient-use system saves sterilization cost and time.



Jazz
Item #1008



Cap
Item #1009



Knit Edge
Item #1010

PROBLEM: Making adjustments to porcelain crowns can be time consuming. Having to send crowns back to the lab for re-glazing is inefficient and requires multiple appointments for the patient.

SOLUTION: Intra-oral polishing allows you to adjust and polish porcelain crowns to their original glaze in just one visit. Save time and improve restoration quality. Jazz™ offers 2 porcelain polishing systems which are also compatible with metals.

When you need to adjust, smooth and polish your porcelain & metal restoration:

JAZZ™ PORCELAIN & METAL
P3S 3-STEP
POLISHING
SYSTEM

Diamond impregnated to achieve brilliant results on porcelain and metal restorations.

A specially formulated rubber matrix allows for an excellent polish with minimal pressure.

Compatible with all porcelain and metal restorations which eliminates the need for multiple polishing systems.

Re-usable for multi-patient use.

1. Use the coarse grit polisher for reducing and trimming.



Coarse Flare
Item #3904



Coarse Cap
Item #3905



Coarse Knife Edge
Item #3906



Coarse Javel Flare
Item #3907

2. Follow with the medium grit polisher for smoothing.



Medium Flare
Item #3908



Medium Cap
Item #3909



Medium Knife Edge
Item #3910



Medium Small Flare
Item #3911

3. Finish with the fine grit polisher to achieve a high luster.



Fine Flare
Item #3912



Fine Cap
Item #3913



Fine Knife Edge
Item #3914



Fine Javel Flare
Item #3915

When you need to save time without compromising quality:

JAZZ™ PORCELAIN & METAL
P2S 2-STEP
SINGLE-PATIENT
POLISHING SYSTEM

Silicon carbide particles embedded in a specially formulated silicone matrix to achieve a high gloss finish.

2-step system smoothes and polishes quickly.

Compatible with all porcelain and metal restorations which eliminates the need for multiple polishing systems.

Single-patient-use system saves sterilization cost and time.

1. Use the medium grit polisher for reducing and smoothing.



Medium Flare
Item #3917



Medium Cap
Item #3918

2. Finish with the fine grit polisher to achieve a high gloss.



Fine Flare
Item #3919



Fine Cap
Item #3920

SWITE™
JAZZ™
POLISHERS

All Jazz™ Polishers are right angle latch-type.

FOR PORCELAIN & METALS

PROBLEM: How can I keep my polishers organized?

SOLUTION: Jazz™ polishing systems are available in autoclavable bar blocks. Our single-patient systems are available in convenient assortment packs.

Kits and Assortment Packs

COMPOSITE POLISHERS

JAZZ
SUPREME

COMPOSITE
3-STEP
UNIVERSAL
POLISHING SYSTEM

4 Pc. Kit in Autoclavable Bar Block Item #80015

Kit includes (3) of each:
#1101 Diamond Flare
#1102 Diamond Cut Edge
#1103 Diamond Cup
#1104 Diamond Seal Flare
Available in Bar Block



JAZZ
C2S

COMPOSITE
2-STEP
POLISHING
SYSTEM

8 Pc. Kit in Autoclavable Bar Block Item #80030

Kit includes (3) of each:
#1105 No-Gloss Flare
#1106 No-Gloss Cut Edge
#1107 No-Gloss Cup
#1108 No-Gloss Seal Flare
#1109 Fluor Flare
#1110 Fluor Cut Edge
#1111 Fluor Cup
#1112 Fluor Seal Flare
Available in Bar Block



JAZZ
C1S

COMPOSITE
1-STEP
SINGLE-PATIENT
POLISHING SYSTEM

12 Pc. Assortment Pack Item #80030

Assortment Pack includes (3) of each:
#1101 Diamond Cup
#1102 Diamond Cut Edge
#1103 Diamond Flare



PORCELAIN & METAL POLISHERS

JAZZ
P3S

PORCELAIN & METAL
3-STEP
POLISHING
SYSTEM

12 Pc. Kit in Autoclavable Bar Block Item #89016

Kit includes (3) of each:
#1113 Crown Flare
#1114 Crown Cup
#1115 Crown Cut Edge
#1116 Crown Seal Flare
#1117 No-Gloss Flare
#1118 No-Gloss Cup
#1119 No-Gloss Cut Edge
#1120 No-Gloss Seal Flare
#1121 Fluor Flare
#1122 Fluor Cup
#1123 Fluor Cut Edge
#1124 Fluor Seal Flare
Available in Bar Block



JAZZ
P2S

PORCELAIN & METAL
3-STEP
SINGLE-PATIENT
POLISHING SYSTEM

12 Pc. Assortment Pack Item #89021

Assortment Pack includes (3) of each:
#1113 No-Gloss Cup
#1114 No-Gloss Flare
#1115 Fluor Cup
#1116 Fluor Flare



WHITE
JAZZ
POLISHERS



1145 Torbin Avenue

Lakewood, NJ 08701, USA

732-902-1100

1-800-632-2877

www.sewhitebus.com

Star® is a trademark of SE White Star Inc.
SE White® is a registered trademark of SE White Star Inc.
All rights reserved.
SE White® is an ISO 15001 registered Company.

SEWE - 01875
2007